



UPPER

Efficient, stable, convenient and rapid operation with full automation

ExoFaster - 600

Exosome Automated Isolation Platform



Efficient stable

Simultaneously isolate high-quality sEVs from 9 samples in 48 minutes



Ultra-high purity

Extracted sEV purity $\geq 1.0 \times 10^8$ particles/ μg protein



Intact morphology

Extracted sEVs show intact saucer-shaped membranes (TEM)



User-friendly

Standardized and user-defined programs

Size-exclusion chromatography

Enables rapid, high-yield isolation of 30–150 nm small extracellular vesicles (sEVs) from biological specimens, including cell culture supernatants and serum/plasma.



SERUM/
PLASMA



URINE



CEREBRO-
SPINAL FLUID



MILK



CELL
SUPERNATANT



PLANT



BACTERIAL
CULTURE


Upper BioTech Pharma Co., Ltd.

London BioScience Innovation Centre, 2 Royal College St, London NW1 0NH
www.upperbiotech.com | general@upperbiotech.com

Function

Model	ExoFaster-600
Function	Standardized and user-defined programs
Sample Throughput	1-9 samples selectable
Purification Duration	Only 48 minutes to purify sEVs from 9 samples
Exosome yield	A single extraction of 1 mL of plasma can yield $\geq 10^{11}$ particles/mL
Exosome proportion	$\geq 90\%$ of particles are in the 40 - 150 nm size range

Performance



Comparison

